

## LGJ-18 Standard Type Experimental Freeze Dryer



LGJ-18 Standard Type Experimental Freeze Dryer is widely used in medicine, pharmacy, biology research, chemical industry and food production, etc. After the freeze drying process, a long term preservation is much easier. They can be restored to the original state and mantain their chemical and biological characteristics after being watered. This model is suitable for experiment in laboratory and small amount production.

## Features:

- . All-in-one structure, small size, no external flange, easy to use, no leakage.
- . All materials in contact with the product use inert materials to meet the requirements of the GLP.
- . The cold trap and console are made of stainless steel, which is anti-corrosion and easy to clean.
- . All stainless steel inflated /water release valve is designed and produced by our company, which is safety, anti- corrosion, no leakage.
- . With sample pre-frozen function, large cold trap opening, no inner coil,and low temperature refrigerator is needless.
- . Patented gas diversion technology, and strong ice holding ability.
- . World- known brand compressors, high efficiency, long life, low noise.
- . Well-known brand vacuum pump with a high pumping speed to achieve a higher ultimate vacuum.
- . Vacuum pump protection function can set the cold trap temperature to protect the vacuum pump life.
- . Professionally designed FD-LAB freeze-drying machine control system + SH-HPSC-I modular controller with high reliability and stability.
- . Intelligent data recording system can real-time record and display the cold trap temperature curve, sample temperature curve, vacuum curve.

## Parameters:

• freeze drying area: 0.18m²/0.27m²

tray size : Φ240mm
No. of tray : 4pcs/6pcs
tray spacing : 70mm

• cold trap temperature : ≤-56C, optional≤-80°C

cold trap dept : 400mmcold trap dia. : Φ270mm

ice collecting capacity : 6kg/24h

pumping speed : 4L/Sultimate vacuum : ≤5pa

power : 1400wweight : 105kg

• main unit dimensions: 630×580×970mm

-80°C main unit dimensions: 810×580×950mm

drying chamber size : Φ300×445mm

cooling mode : wind cooleddefrosting mode : air cooled

material capacity: 1.8L/2.7L (thickness 10mm)

Optional accessories and functions

cold trap temperature -80°C

vacuum pump anti-returned oil device

• ·cold trap electrical heat defrosting

• ·oil mist filter system

• ·imported vacuum pump

non-standard customized